RESEARCH

BMC Oral Health

Open Access

Color stability of two different resin matrix ceramics: randomized clinical trial



Aliaa Ibrahim Mahrous^{1,2*}, Aya A. Salama², Alshaimaa Ahmed Shabaan³, Ahmed Abdou^{4,5}, and Mohamed Mostafa Radwan⁶

Abstract

Background One of the most common causes of aesthetic failure and restoration replacement is the tooth restorations color mismatch specifically after aging.

Methods One hundred and two participants with endodontically treated first molar were selected clinically. The patients were randomly splited into two groups and restored either with Cerasmart hybrid ceramic or Vita Enamic polymer infiltrated ceramic network (PICN) crowns and cemented using dual cure adhesive cement. The color difference (ΔE) values after cementation at 0 (Baseline), 6, and 12 months of use were obtained by quantification of L*, a*, and b* values with a digital spectrophotometer. Mann–Whitney test used to compare between tested groups at each time point and between (α =0.05).

Results At 6 months follow-up intervals, Vita Enamic group showed the highest significant ΔL^* (p = 0.035) and Δa^* (p < 0.001) compared to Cerasmart group. Δb^* and ΔE showed no significant difference between both groups (p > 0.05). Furthermore, all color parameters of both groups showed statistically significant difference at 12 months follow-up intervals. After 12 months, Vita Enamic restorations presented higher color change compared to Cerasmart restorations with a (p-value of 0.0120). When comparing the total color difference ΔEt through-out the follow-up intervals of Vita Enamic & Cerasmart groups, there were insignificant difference (p = 0.263).

Conclusion Both hybrid materials demonstrated comparable color stability after 1 year of clinical service within clinical acceptance range. However, Cerasmart demonstrated a better colour stability after 1 year.

Trial registration ClinicalTrials.gov (ID: NCT05501808) 15/8/ 2022- 'retrospectively registered'.

Keywords Dental restoration, Ceramics, Resin matrix, Color stability, Spectrophotometer

*Correspondence:

- Aliaa Ibrahim Mahrous
- aim11@fayoum.edu.eg
- ¹ Fixed Prosthodontic Department, Faculty of Dentistry, Fayoum University, Fayoum, Egypt
- ² Fixed Prosthodontics Department, Faculty of Dentistry, October
- University for Modern Sciences and Art, 6th of October City, Egypt
- ³ Faculty of Dentistry, Fayoum University, Fayoum, Egypt
- ⁴ Prosthetic Dentistry Department, Biomaterials Division, Faculty of Dentistry, King Salman International University, El Tur, South Sinai, Egypt

⁵ Faculty of Dentistry, Al-Ayen University, Thi-Qar, Iraq

⁶ Fixed Prosthodontic Department, Faculty of Dentistry, Bani Suif University, Bani Suif, Egypt

Background

The aesthetic CAD/CAM materials have undergone significant advancements, changing the treatment from a two-step, bi-layered, high-strength ceramic restoration to a single step, monolithic restoration that eliminates veneering chipping and fracture issues [1, 2]. Obtaining the ideal optical qualities of natural teeth with artificial restorative materials is one of modern dentistry's primary challenges. Nowadays, fixed dental prosthesis can be made from a wide range of dental materials. Ceramics and modified composites are largely in question for this purpose, as they



© The Author(s) 2023. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.gr/licenses/by/4.0/. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated in a credit line to the data.

Content courtesy of Springer Nature, terms of use apply. Rights reserved.

are manufactured to closely resemble natural teeth in color, which is a desirable trait for many patients [3–6]. Ceramics have been the preferred material due to their exceptional aesthetics, high strength, and biocompatibility. However, ceramic restorations are challenging as it can cause opposing teeth to wear out excessively and have poor repair options [6].

With the advancements in adhesive dentistry and resin composites, a new unique material that has been referred to as "nano-ceramics" or "hybrid ceramics" in the commercial market was developed in the last ten years with high chance of success as aesthetic restoration [7]. A shock absorber hybrid ceramic is a modern alternative to traditional ceramics for use in dental restorations such as inlays, partial crowns, and full crowns [8, 9] with physical qualities highly comparable to those of teeth [10]. Hybrid ceramics have properties between those of ceramics and composites [11, 12], including somewhat high fracture toughness, elastic modulus, hardness, and rigidity. Hybrid ceramics are more flexible, less expensive, and less prone to fracture. Also, when compared to all ceramic materials and composites, it gives rise to reduced abrasion of opposing natural teeth with an important advantage over ceramics [13-17]. In addition to easier fabrication, being repaired easily and used right in the dentist's chair without needing to be fired first [11, 18].

Hybrids are considered as infiltrating porous ceramic with polymers. There are two categories on the market based on microstructure; dispersed filler (DF) and polymer-infiltrated ceramic network (PICN). The ceramic and polymer phases of PICN, also known as resin interpenetrating network ceramic, form a continuous three-dimensional interconnected network [8]. A PICN example is Vita Enamic (VITA Zahnfabrik) while GC America's Cerasmart is an example for resin with dispersed fillers [6, 9, 19].

The first material in this classification to be introduced is Vita Enamic. The product's attributes, according to the manufacturer, lessen the likelihood of crack propagation and fracture by combining the advantages of composites and ceramics to relay strength and elasticity and behave like dentin. Its flexural strength ranges from 150 to 160 MPa and its composition is 86% ceramic and 14% polymer. Due to the special ceramic composition, it is unique in that it can be etched with 5% hydrofluoric acid for 60 s on the intaglio surface [8, 9].

Cerasmart nanoceramic block was introduced to combine the advantageous qualities of ceramic and resin technology. The material is made up of 29% resin and 71% silica and barium glass nanoparticles. This material, which has a flexural strength of 238 MPa, is intended for Page 2 of 8

anterior, posterior, and implant restorations that require the least amount of preparation [6, 10].

Aesthetic considerations in dental rehabilitation highlight the importance of the color stability of dental materials. Different foods and beverages can stain dental restorations when they come into contact. All these may cause the restorations to fail aesthetically. Composites have a far higher discoloration potential when compared to ceramics [20, 21].

The amount of discoloration is mostly determined by the composition of the material, type of environment, and contact period [22]. The lightness, Hue, and chroma components of the traditional CIELab color formula are denoted by the letters L*, H*, and C respectively. Accordingly, L*, a* and b* where L* stands for lightness and takes a number [11] between 0 to 100, with 0 refers to black while 100 stands for white. The saturation on the letter a* is the saturation on the redto-green axis, and b* blue-yellow axis. In order to calculate color differences, a formula is used to express these parameters as "E values" [23, 24].

In clinical situations with different dental materials, the evaluation of color distinctions is quite important. The color stability of resin matrix ceramics is critical, as clinically observable over time. Due to the high demand for cosmetic procedures, clinicians must be careful when selecting restorative materials, since this is one of the most important elements determining long-term treatment effectiveness [25].

The durability of color in contemporary hybrid ceramics is mostly unknown. The current in-vivo study aimed to probe the color stability of two hybrid ceramics (Vita Enamic; EN) and (Cerasmart; CS) under the variety of clinical conditions. The null hypothesis of the present study was that there will be no statistically significant difference regarding color stability of the tested two resin matrix ceramics.

Materials & methods

The study was carried out at the Fixed Prosthodontics Department, Faculty of Dentistry, Fayoum University, between June 2021 and August 2022. The Declaration of Helsinki, which governs research involving human subjects, was followed [26]. The Fayoum University's supreme ethics committee EC2140 gave the study protocol a thorough review, approval, and retrospective registration in the ClinicalTrials.gov database. (Identifier: NCT05501808) on 15th of August 2022. Prior to the study's start, each participant provided their written, informed consent to share their data, and they were free to leave the study at any time without having to give a reason.

Sample size

Based on previous research [27], a sample size calculation was made using STATA16 with a desired alpha of 0.05 and power of 80% to distinguish between the null hypothesis, which states that the means of both groups are 1.722, and the alternative hypothesis, which states that the means of the Vita Enamic hybrid ceramic group are 2.2, with estimated group standard deviations of 0.854 and the ratio of the sample sizes in groups to be 1. It was determined that 90 participants would make up the entire sample for this study. However, 102 participants (51 in each group), presuming potential dropouts, were included in this study.

Study design and randomization

To ensure balance in the number of participants assigned to each group, this study was prepared as a doubleblinded prospective randomized controlled trial. The randomization was unstratified and carried out via a random block design with blocks of size 2, 4, and 6 to ensure balance in the number of patients assigned to each group. The study was designed in accordance with the Consolidated Standards of Reporting Trials (CONSORT) statement of 2012 (Fig. 1) [28].

One hundred two participants with previous endodontically treated first molar were eligible for full coverage hybrid ceramic crown in lower 1st molar. The inclusion criteria were as follows: (1) Age of the participants from 18–60 years old who can read and sign informed consent document. (2) Participants with proper endodontic treatment (3) participants with no systemic conditions that could affect the procedures. (4) No active periodontal disease. (5) Participants with stable normal occlusion. While the exclusion Criteria were; (1) Participants with high caries index. (2) Participants with smoking habits. (3) Participants with parafunctional habits as clenching/bruxism. (4) Participants with drugs that affect oral health. (5) Pregnant females.

The participants were allocated to one of the two groups based on hybrid ceramic type used, (Table 1). Vita Enamic group (n=51): crowns were milled from (Vita Enamic, Vita Zahnfabrik, Bad Säckingen, Germany)




Fig. 1 Consort flow chart

Table 1 Desc	cription of h	nybrid	ceramics	utilized	in the study	/

Hybrid system	Composition	Manufacturer	
GC Cerasmart (Resin matrix ceramic)	Resin matrix: Bis-MEPP, UDMA, dimethacrylate	GC Corp., Tokyo, Japan	
	Inorganic filler: silica, barium glass (71.0 wt%)		
Vita Enamic (Resin matrix ceramic (polymerinfiltrated-ceramic-network (PICN))	Resin matrix: UDMA, TEGDMA	Vita-Zahnfabrik, Bad	
	Inorganic filler: feldspar ceramic enriched with aluminum oxide (86.0 wt%)	Sackingen, Germany	

TEGDEMA Triethylene glycol dimethacrylate, Bis-MEPP (2,2-Bis (4-methacryloxypolyethoxyphenyl) propane, UDMA urethane dimethacrylate)

while Cerasmart group (n=51): crowns were milled from (Cerasmart; CS,GC Corp., Tokyo, Japan). The material used served as the primary predictor. The assessors and participants were both blinded for the duration of the study.

Intervention

All participant underwent clinical, radiographic evaluation of all cases with shade selection using Vita System 3-D Master (VITA) was performed. During the preparatory visit, any carious lesions were removed and restored using a core build-up restorative material (GC Fuji lining LC; GC Corp, Tokyo, Japan) upon requirements. All teeth were prepared in a standardized reduction parameters following the manufacturer's recommendations with a deep chamfer finish line, 1–1.5 mm occlusal reduction and 0.8–1.5 mm axial reduction. Astringent retraction paste (VOCO GmbH) was used into sulcus prior to digital impression.

CEREC, Dentsply Sirona was used to scan and construct the entire crowns in current study. An intra-oral impression was captured using a powder-free intraoral scanner CEREC Omnicam (Dentsply Sirona) following the manufacturer's recommendations. The prepared tooth, opposing teeth and interocclusal relationship was recorded and the biogeneric individual design mode of CAD software was used for the design. CEREC MCXL milling machine (Dentsply Sirona) was used for CAM fabrication of both groups. Checking of all crowns for proper fit with adjustments to occlusal and proximal contacts were made followed by final finishing with NTI CeraGlaze (Axis, Coppell) on an electric handpiece while following the manufacturer's instructions at 10,000 rpm for 10 s with a fine finishing wheel and at 5,000 rpm for 10 s with a high-shine polishing wheel.

At delivery visit, all prepared teeth followed same protocol for cementation. All prepared teeth were cleaned using Pumice and Robinson brush before isolated with rubber dam. Phosphoric acid etchant 37% was applied for 20 s, washed and air dried & bonding agent was applied. The fitting surfaces of all crowns were etched 60 s using 5% hydrofluoric acid according to the manufacturer's instructions. The etched surfaces were thoroughly water sprayed and oil-free dried before application of silane coupling agent for 60 s. Dual-cure resin cement (Vriolink Esthetic, Ivoclar Viadent AG, Schaan, Liechtenstein) was used. Using a sharp explorer, extra cement was scraped off. Air-blocking gel (Oxiguard II: Kuraray Noritake Dental Inc.: Tokyo, Japan) was finally applied and light-cured (Elipar Deepcure-S: 3 M Espe: St. Paul, USA) about 20 s.

A spectrophotometer (VITA Easy-shade Advance 4.0, VITA Zahnfabrik, Bad Säckingen, Germany) was used to assess the color of all participants at the follow-up appointments (baseline = 24 h after cementation, 6 months and 12 months). It was calibrated before each measurement. CIELAB parameters (L*, a*, b*) were recorded in at each follow-up appointment with the same operator who was blinded. L* indicated the material lightness, a* denotes hues change across the red-green axis, while b* showed the hues difference across the yellow-blue axis. Finally, ΔL^* , Δa^* , Δb^* and ΔE were calculated for 6, 12 month and finally the overall color changes over the year. CIELAB (ΔE_{ab}) was used for calculating the color differences as follows;

$$\begin{split} \Delta E_6 &= [(\Delta L_6 = L_6^* - L_0^*)^2 + (\Delta a_6 = a_6^* - a_0^*)^2 + (\Delta b_6 = b_6^* - b_0^*)^2]^{1/2} \\ \Delta E_{12} &= [(\Delta L_{12} = L_{12}^* - L_6^*)^2 + (\Delta a_{12} = a_{12}^* - a_6^*)^2 + (\Delta b_{12} = b_{12}^* - b_6^*)^2]^{1/2} \\ \Delta E_t &= [(\Delta L_t = L_{12}^* - L_0^*)^2 + (\Delta a_t = a_{12}^* - a_0^*)^2 + (\Delta b_t = b_{12}^* - b_0^*)^2]^{1/2} \end{split}$$

Statistical analysis

Data checked for normality using Kolmogorov Smirnova test. Data displayed non-normal distribution, so Mann–Whitney test used to compare between tested group at each time point and between time points for each group. Significant level was set at p=0.05. Statistical analysis was done using IBM SPSS (version 26, Armonk, NY, USA).

Results

This study included 102 participants allocated into 2 groups equally. No complications were observed and no drop out except for one participant in each group which was excluded from the analysis. At 6 months follow-up intervals, Vita Enamic group showed the highest significant ΔL^* , (p=0.035) and Δa^* , (p<0.001) compared to Cerasmart group. While, Δb^* , and ΔE showed insignificant difference between both groups (p>0.05). Furthermore, all color parameters of Cerasmart showed lower significant values compared to Vita Enamic at 12 months follow-up intervals.

Regardless of the restoration materials, all color parameters showed significant decrease except for ΔL^* , within Cerasmart group (Table 2 and Fig. 2).

When comparing the total color difference ΔE_t through-out the follow-up intervals of Vita Enamic & Cerasmart groups, there were insignificant difference (*p*=0.263) (Fig. 3).

Discussion

Over the past few years, there has been a growing trend of patients requesting or expecting high aesthetic restoration. Color stability is an essential clinical aspect when it comes to the aesthetics which in turns affect the

Table 2 Mean [95% CI] for color parameters (Δ L, Δ a, Δ b, and Δ E) at follow-up intervals for different tested groups

Color parameter	Material	6 months	12 Months	p-value	
Δ.	EN	0.3[0.3 to 0.8]	0.1[0.0 to 0.3]	< 0.001	
	CS	0.1[-1.2 to 1.1]	0.0[0.0 to 0.3]	0.8488	
p-value		0.0351	0.0216		
Δα	EN	0.3[0.1 to 0.8]	0.0[-0.1 to 0.1]	< 0.001	
	CS	-1.1[-1.1 to -0.3]	0.0[0.0 to 0.0]	< 0.001	
p-value		< 0.001	0.0044		
Δb	EN	0.3[0 to 0.5]	0.0[-0.1 to 0.1]	< 0.001	
	CS	0.3[0.3 to 0.5]	0.0[0.0 to 0.0]	< 0.001	
p-value		0.1896	0.0044		
ΔΕ	EN	1.2[0.7 to 2]	0.3[0.1 to 0.8]	< 0.001	
	CS	1.7[1.2 to 2.1]	0.0[0.0 to 1.2]	< 0.001	
p-value		0.0853	0.0120		

Significant level was set at p < 0.05




Fig. 2 Box Plot showing comparison between Vita Enamic and Cerasmart in all Color parameters (ΔL^* , Δa^* , $\Delta b^* \& \Delta E$) at follow-up intervals of 6month and 12 month

long-term clinical success of resin-based ceramic restorations. Plaque accumulation, solution stains, surface roughness, and chemical reactions can all affect the color of these materials [29, 30].

As instrumental color measurement has the advantages of being objective and quantified, CIE lab formula was selected in the current study to determine any color changes [31]. The color change values were represented by ΔE which is defined as the numerical distance between the L*a*b* coordinates of two colors. It




Fig. 3 Box Plot showing the ΔE_t (total color difference after 1 year)

was reported that there were a range before detectable color changes. $\Delta E < 1$ were undetectable by the human eye; $1.0 < \Delta E < 3.3$ was seen only by a skilled person and considered clinically acceptable while, $\Delta E > 3.3$ was considered easily observed and not acceptable clinically [31, 32].

Previous research on hybrid materials has mostly concentrated on mechanical qualities [33], making the characterization of optical properties rather limited. Based on the assumption that resin-based restorative materials are predicted to be less color-stable upon aging due to polymer matrix decomposition, unreacted monomers of the polymerization agents, and extrinsic staining agents, materials with comparable composition and fundamental microstructures were chosen for the current study [5, 34, 35]. Cerasmart and Vita Enamic are both categorised as hybrid materials with different compositions which justifies the necessity to study variations in color stability. Cerasmart material is made up of 29% crosslinked polymer blends and 71% ceramic-like inorganic silicate glass fillers. Contrarily, Vita Enamic is containing 86% feldspathic ceramic and 14% cross-linked polymer blends [36].

Currently, there is no published clinical trial research on the color stability of resin-ceramic infiltrated materials. Therefore, this study compared the color stability of two resin hybrid ceramic restorations in randomized clinical trial. Founded on the study outcomes, the null hypothesis stated was accepted since the color stability of Vita Enamic was comparable to that of Cerasmart hybrid ceramic; with no statistically significant difference after one year of clinical service and all was within the clinically acceptable range of color difference ($\Delta E > 3.3$) [21, 32]. In contrast, Barutçug et al. [37] found that a discoloration greater than the clinically acceptable threshold level ($\Delta E00 = 2.25$) was seen for both Cerasmart and Vita Enamic after a month of exposure to beverages.

The results of the current study compared color alterations in two dental materials over a period of 6, 12 months and the total change. The findings showed that the Vita Enamic group significantly higher than the Cerasmart group in terms of ΔL^* (lightness) and, Δa^* (hue) values at 6-month intervals, whereas there was no significant difference in Δb^* (color tone) and E (overall color difference) values between both groups with a *p*-value of 0.0853. This suggested that Vita Enamic in the short term have a more pronounced impact on the teeth's brightness (lightness) and color (hue) than Cerasmart.

Difference in materials color stability across the research timeline could be endorsed to differences in composition and microstructure [36]. The higher color change values of Vita Enamic could be credited to its alumina content (8.31 wt%), that could cause decrease in

translucency with less color stability compared to Cerasmart on the long term of use [38]. Contrary, Cerasmart has a lower proportion of alumina and increased ratio of zirconia and silica nanoparticles surrounded in a highly cross-linked resin matrix, which may explain its greater translucency and color stability [39]. Furthermore, the ageing process disrupts the chemical link between the filler and the resin matrix, which contributes to color instability [40].

Additionally, the sustainability of a restoration's color is mainly influenced by the type of resin matrix, as water absorption by the resin component of the material is a major contributing factor to color changes [7, 41]. Hydrophilic bisphenol A-glycidyl methacrylate (BISGMA) is a methacrylate monomer with high strength and durability that is used commonly [42, 43]. It can cause discoloration due to the high susceptibility to leaching from the material. Paolone et al. [44] suggested also that this monomer is responsible for water uptake and, therefore, possibly for discoloration. Hydrophobic urethane dimethacrylate (UDMA) is another methacrylate monomer that has high strength and low shrinkage. It may explain the color stability of hybrid ceramics by decreasing the polymerization shrinkage and stress so reducing color changes over time. TEGDMA is known for its low viscosity and ability to decrease the curing time but it has been shown to be more prone to leaching out of dental composites compared to other monomers, which can potentially cause discoloration or color shift over time [42–45].

Saba et al. [32] mentioned that increasing the TEG-DMA proportion from 0 to 1% will increase the water uptake of Bis- GMA based resins from 3 to 6% respectively. The high wt% of TEGDMA in Vita Enamic will result in increased water sorption and thus allow the penetration of any hydrophilic colorant into the resin matrix [40, 43]. Although, UDMA is hydrophobic compared to BIS-GM, di-methacrylates, can create cross-linked networks that trap unreacted monomers with plasticizers and the formation of a more porous structure. This can increase the water sorption of the material [32]. This may explain the more color changes of Vita Enamic compared to Cerasmart (*p*-value of 0.0120) especially after 12 month of oral environment aging of both materials.

In the same context, Sarikaya et al. [46] found that Vita Enamic had a significantly greater color difference than other resin nano ceramics. Stamenković et al. [47] noticed no differences in the color stability of different materials. Oppositely, Arif et al. [42] mentioned that nanoceramic materials can be best avoided for veneers in areas with high aesthetic requirements especially by coffee drinkers. Also, Al Amri et al. [43] compared color stability of five CAD/CAM materials (Lava Ultimate (LU), Cerasmart (CS), Vita Enamic (VE), Crystal Ultra (CU) IPS e.max (IPS) after soaking in coffee drink resulted in observable color changes except for PICN materials which revealed acceptable color changes but, resin nanoceramic samples has the biggest color differences, with Cerasmat showing greater color changes (Δ E00=2.09) than Vita Enamic (Δ E00<1.8).

In addition, Materials can experience color changes due to many factors such as abrasion, erosion, exposure to saliva, food, and drinks [44, 48] The ability of a material to withstand occlusal wear is a crucial aspect for the long-term success of a prosthesis. In a study analyzing the wear pattern of resin ceramic matrix, Vita Enamic demonstrated a higher wear pattern compared to Cerasmart. The filler system, including both the type and size, has undergone significant improvements in the mechanical and physical properties of resin matrix materials [16]. Smaller filler sizes resulted in better wear behavior, while larger filler sizes led to a more protruded surface with higher risk of scratching and subsequently more color changes [12, 17]. The manufacturing procedures used to create Vita Enamic and Cerasmart could also have an impact on how stable their colors are. For instance, the structure and properties of the materials may be impacted by the curing conditions, temperature, and pressure used during the manufacturing process [49].

At the end, both groups revealed statistically insignificant variations in color stability, proving that both materials underwent color change. This might imply that both materials' color stability would be compromised over a longer time frame. Even though the differences may be slight and invisible to patients or clinicians and within clinically acceptable range [10, 11].

Limitations

The first drawback of this study is the small sample size, which can be rationalized by the fact that it was calculated using the 5% margin of error assumption. Future studies should take lower-margin errors into consideration. Other limitations may be that clinical color change may be impacted by routine toothbrushing with various kinds of toothpaste. It is advised that more research be done to assess the optical characteristics of hybrid materials at various thicknesses and under several aging procedures. A longer-term clinical investigation should confirm these clinical findings.

Conclusions

- 1. Cersmart demonstrated better color stability after 1 year of clinical service.
- 2. Both Vita Enamic and Cerasmart undergo color changes over time within clinical acceptance range.

Abbreviations

- µ Mean
- SD Standard deviation
- CI Confidence level
- Sig. Significance level

Acknowledgements

Not applicable.

Authors' contributions

AM: study design, conceptualization and major contributor in writing the manuscript- MR: conceptualization and interpreted the patient data- AAS: interpretation of data- AA: data analysis and article revision- AS: conceptualization, interpreted the patient data and drafting of the article.

Funding

Open access funding provided by The Science, Technology & Innovation Funding Authority (STDF) in cooperation with The Egyptian Knowledge Bank (EKB). This research received no specific grant from any funding agency in the public, commercial, or non-profit sectors.

Availability of data and materials

The datasets used and/or analysed during the current study are available from the corresponding author on reasonable request. For privacy reasons, however, individual data allowing for the identification of participants cannot be made available.

Declarations

Ethics approval and consent to participate

The study was conducted according to the guidelines of the Declaration of Helsinki and approved by the Ethics Committee of Fayoum University. All involved subjects gave their written informed consent before study inclusion. All involved Participants gave their written informed consent for the publication.

Consent for publication

Not applicable.

Competing interests

The authors declare that they have no competing interests.

Received: 15 May 2023 Accepted: 27 August 2023 Published online: 14 September 2023

References

- 1. Suárez MJ, Lozano JF, Paz Salido M, et al. In vitro marginal fit of CAD-CAM monolithic ceramic restorations. BMC Oral Health. 2021;21(1):256.
- Özcan M, Kurbad A, Gurgan S, et al. Clinical outcomes of lithium disilicate versus zirconia single crowns prepared by digital workflows: a multicenter RCT with 3-year follow-up. BMC Oral Health. 2021;21(1):347.
- Higa M, Kano T, Oda Y, et al. Color changes of monolithic zirconia and lithium disilicate after accelerated aging. BMC Oral Health. 2021;21(1):117.
- Wang T, Cao Z, Li Q, et al. Comparison of the color stability of different CAD/CAM materials after immersion in coffee. BMC Oral Health. 2021;21(1):115.
- Acar O, Yilmaz B, Altintas SH, Chandrasekaran I, Johnston WM. Color stainability of CAD/CAM and nanocomposite resin materials. J Prosthet Dent. 2016;115:71–5.
- Abdou A, Takagaki T, Alghamdi A, Tichy A, Nikaido T, Tagami J. Bonding performance of dispersed filler resin composite CAD/CAM blocks with different surface treatment protocols. Dent Mater J. 2021;40(1):209–19.
- Zhang Y, Wang Y, Li Q. Effects of different surface treatments on the color stability of monolithic zirconia ceramics. BMC Oral Health. 2021;21(1):306.
- Coldea A, Swain MV, Thiel N. Mechanical properties of polymer-infiltrated-ceramic-network materials. Dent Mater. 2013;29:419–26.

- He LH, Swain M. A novel polymer infiltrated ceramic dental material. Dent Mater. 2011;27:527–34.
- 10. Seyidaliyeva A, Rues S, Evagorou Z, Hassel AJ, Rammelsberg P, Zenthöfer A. Color stability of polymer-infiltrated-ceramics compared with lithium disilicate ceramics and composite. J Esthet Restor Dent. 2020;32:43–50.
- 11. Yerliyurt K, Sarıkaya I. Color stability of hybrid ceramics exposed to beverages in different combinations. BMC Oral Health. 2022;22(1):180.
- Radwan MM, Abdou A, Nagi N, Younes D, Mahrous A. Wear characteristics of esthetic resin matrix and zirconia reinforced lithium disilicate CAD/ CAM materials: two-body wear and surface topography analysis. Minerva Dent Oral Sc. 2023. https://doi.org/10.23736/S2724-6329.23.04779-4.
- Ashtiani AH, Azizian M, Rohani A. Comparison the degree of enamel wear behavior opposed to polymer-infiltrated ceramic and Feldspathic porcelain. Dent Res J. 2019;16:71–5.
- 14. Yin R, Kim YK, Jang YS, et al. Comparative evaluation of the mechanical properties of CAD/CAM dental blocks. Odontology. 2019;107:360–7.
- Choi JW, Song EJ, Shin JH, Jeong TS, Huh JB. In vitro investigation of wear of CAD/CAM polymeric materials against primary teeth. Materials. 2017;10:1410–22.
- Xu Z, Yu P, Arola DD, Min J, Gao S. A comparative study on the wear behavior of a Polymer Infiltrated Ceramic Network (PICN) material and tooth enamel. Dent Mater. 2017;33:1351–61.
- Zhi L, Bortolotto T, Krejci I. Comparative in vitro wear resistance of CAD/CAM composite resin and ceramic materials. J Prosthet Dent. 2016;115:199–202.
- Alomari QD, Alkhraisat MH, Al-Wahadni AM, Al-Sabbagh R, Al-Hiyasat AS, Al-Jamal HA. Fracture resistance and mode of failure of CAD/ CAM resin composites with different translucencies. BMC Oral Health. 2021;21(1):195.
- Skorulska A, Piszko P, Rybak Z, Szymonowicz M, Dobrzyński M. Review on polymer, ceramic and composite materials for CAD/CAM indirect restorations in dentistry-application, mechanical characteristics and comparison. Materials (Basel). 2021;14(7):1592.
- 20. Samra AP, Pereira SK, Delgado LC, et al. Color stability evaluation of esthetic restorative materials. Braz Oral Res. 2008;22:205–10.
- Faris TM, Abdulrahim RH, Mahmood MA, Mhammed Dalloo GA, Gul SS. In vitro evaluation of dental color stability using various aesthetic restorative materials after immersion in different drinks. BMC Oral Health. 2023;23(1):49.
- 22. Alsilani RS, Sherif RM, Elkhodary NA. Evaluation of colour stability and surface roughness of three CAD/CAM materials (IPS e.max, Vita Enamic, and PEEK) after immersion in two beverage solutions: an in vitro study. Int J Appl Dent Sci. 2022;8:439–49.
- 23. International Commission on Illumination. CIE015:2004: Technical Report: Colorimetry. Vienna: 2004.
- 24. Commision Internationale de l'Eclairage (CIE). http://cie.co.at/publicatio ns/colorimetry-part-6-ciede2000-colour-diference-formula. Accessed 03 Feb 2021.
- Gugelmin BP, Miguel LCM, Baratto Filho F, Cunha LFD, Correr GM, Gonzaga CC. Color stability of ceramic veneers luted with resin cements and pre-heated composites: 12 months follow-Up. Braz Dent J. 2020;31(1):69–77.
- 26. Mundial AM. Declaración de Helsinki de la AMM Principios éticos para las investigaciones médicas en seres humanos. 2019.
- 27. Kanat-Ertürk B. Color stability of CAD/CAM ceramics prepared with different surface finishing procedures. J Prosthodont. 2020;29(2):166–72.
- Schulz KF, Altman DG, Moher D. CONSORT 2010 statement: updated guidelines for reporting parallel group randomised trials. Trials. 2010;11(1):1–8.
- Alqahtani MQ, Chinelatti MA. Color stability of CAD/CAM materials versus direct resin composites: a systematic review. BMC Oral Health. 2019;19(1):5.
- Mohamed I, Ehab F, Rasha B. A comparison of color stability between hybrid ceramic and veneers: an in vitro study. Int Arab J Dent. 2019;10(1):3.
- Jiang T, Ma L, Lv P, Zhao Y, Zhang Y. The effect of different staining solutions on the color stability of dental zirconia ceramics before and after polishing. J Prosthet Dent. 2021;126(4):576–83.
- Saba DA, Salama RA, Haridy R. Effect of different beverages on the color stability and microhardness of CAD/CAM hybrid versus feldspathic ceramic blocks: an in-vitro study. Futur Dent J. 2017;3:61–6.

- Albero A, Pascual A, Camps I, et al. Comparative characterization of a novel CAD/CAM polymer-infiltrated-ceramic-network. J Clin Exp Dent. 2015;7:e495–500.
- Awad D, Stawarczyk B, Liebermann A, et al. Translucency of esthetic dental restorative CAD/CAM materials and composite resins with respect to thickness and surface roughness. J Prosthet Dent. 2015;113:534–40.
- Facenda JC, Borba M, Corazza PH. A literature review on the new polymer-infiltrated ceramic-network material (PICN). J Esthet Restor Dent. 2018;30:281–6.
- 36. Gawriolek M, Sikorska E, Ferreira LF, et al. Color and luminescence stability of selected dental materials in vitro. J Prosthodont. 2012;21:112–22.
- Barutçugil Ç, Bilgili D, Barutcigil K, Dündar A, Büyükkaplan US, Yilmaz B. Discoloration and translucency changes of CAD-CAM materials after exposure to beverages. J Prosthet Dent. 2019;122(4):325–31.
- Sen N, Us YO. Mechanical and optical properties of monolithic CAD-CAM restorative materials. J Prosthet Dent. 2018;119(4):593–9.
- 39. Lee YK. Influence of scattering/absorption characteristics on the color of resin composites. Dent Mater. 2007;23:124.
- Ferracane JL. Hygroscopic and hydrolytic effects in dental polymer networks. Dent Mater. 2006;22:211–22.
- Liu C, Zhang Y, Wang Y. Effects of resin matrix on translucency and color stability of dental ceramics. J Prosthodontic Res. 2021;65(2):185–91.
- 42. Arif R, Yilmaz B, Johnston WM. In vitro color stainability and relative translucency of CAD-CAM restorative materials used for laminate veneers and complete crowns. J Prosthet Dent. 2019;122:160–6.
- Al Amri MD, Labban N, Alhijji S, Alamri H, Iskandar M, Platt JA. In vitro evaluation of translucency and color stability of CAD/CAM polymer-infiltrated ceramic materials after accelerated aging. J Prosthodont. 2021;30:318–28.
- Paolone G, Mandurino M, De Palma F, Mazzitelli C, Scotti N, Breschi L, Gherlone E, Cantatore G, Vichi A. Color stability of polymer-based composite CAD/CAM blocks: a systematic review. Polymers (Basel). 2023;15(2):464.
- Peutzfeldt A, Garcia-Godoy F, Flury S, Eliades G, Lussi A. Influence of coating procedure on the performance of resin composite restorations after 2 years. Clin Oral Investig. 2019;23(8):3165–317.
- Sarıkaya I, Yerliyurt K, Hayran Y. Effect of surface finishing on the colour stability and translucency of dental ceramics. BMC Oral Health. 2018;18:40.
- Stamenković DD, Tango RN, Todorović A, Karasan D, Sailer I, Paravina RD. Staining and aging-dependent changes in color of CAD-CAM materials. J Prosthet Dent. 2021;126(5):672–8.
- Mansour CH, Kheiralla LS, Ghanem LA. Color stability and surface roughness of CAD/CAM ceramill composite and cerasmart endocrowns. J Dent Med Sci. 2018;17:13.
- Kocak E, Ekren O, Johnston W, Ucar Y. Analysis of color differences in stained contemporary esthetic dental materials. J Prosthet Dent. 2021;126:438–45.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Terms and Conditions

Springer Nature journal content, brought to you courtesy of Springer Nature Customer Service Center GmbH ("Springer Nature").

Springer Nature supports a reasonable amount of sharing of research papers by authors, subscribers and authorised users ("Users"), for smallscale personal, non-commercial use provided that all copyright, trade and service marks and other proprietary notices are maintained. By accessing, sharing, receiving or otherwise using the Springer Nature journal content you agree to these terms of use ("Terms"). For these purposes, Springer Nature considers academic use (by researchers and students) to be non-commercial.

These Terms are supplementary and will apply in addition to any applicable website terms and conditions, a relevant site licence or a personal subscription. These Terms will prevail over any conflict or ambiguity with regards to the relevant terms, a site licence or a personal subscription (to the extent of the conflict or ambiguity only). For Creative Commons-licensed articles, the terms of the Creative Commons license used will apply.

We collect and use personal data to provide access to the Springer Nature journal content. We may also use these personal data internally within ResearchGate and Springer Nature and as agreed share it, in an anonymised way, for purposes of tracking, analysis and reporting. We will not otherwise disclose your personal data outside the ResearchGate or the Springer Nature group of companies unless we have your permission as detailed in the Privacy Policy.

While Users may use the Springer Nature journal content for small scale, personal non-commercial use, it is important to note that Users may not:

- 1. use such content for the purpose of providing other users with access on a regular or large scale basis or as a means to circumvent access control;
- 2. use such content where to do so would be considered a criminal or statutory offence in any jurisdiction, or gives rise to civil liability, or is otherwise unlawful;
- 3. falsely or misleadingly imply or suggest endorsement, approval, sponsorship, or association unless explicitly agreed to by Springer Nature in writing;
- 4. use bots or other automated methods to access the content or redirect messages
- 5. override any security feature or exclusionary protocol; or
- 6. share the content in order to create substitute for Springer Nature products or services or a systematic database of Springer Nature journal content.

In line with the restriction against commercial use, Springer Nature does not permit the creation of a product or service that creates revenue, royalties, rent or income from our content or its inclusion as part of a paid for service or for other commercial gain. Springer Nature journal content cannot be used for inter-library loans and librarians may not upload Springer Nature journal content on a large scale into their, or any other, institutional repository.

These terms of use are reviewed regularly and may be amended at any time. Springer Nature is not obligated to publish any information or content on this website and may remove it or features or functionality at our sole discretion, at any time with or without notice. Springer Nature may revoke this licence to you at any time and remove access to any copies of the Springer Nature journal content which have been saved.

To the fullest extent permitted by law, Springer Nature makes no warranties, representations or guarantees to Users, either express or implied with respect to the Springer nature journal content and all parties disclaim and waive any implied warranties or warranties imposed by law, including merchantability or fitness for any particular purpose.

Please note that these rights do not automatically extend to content, data or other material published by Springer Nature that may be licensed from third parties.

If you would like to use or distribute our Springer Nature journal content to a wider audience or on a regular basis or in any other manner not expressly permitted by these Terms, please contact Springer Nature at

onlineservice@springernature.com